DCM045-192-4  
Ed. 01/2015

## T4 ELISA

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## BIBLIOGRAFIA

1. Barker, S.B., "Determination of Protein Bound Iodine" Journal Biological Chemistry, 173, 175 (1948).
2. Chopra, I.J., Solomon, D.H., and Ho, R.S., "A Radioimmunoassay of Thyroxine, "J.Clinical Endocrinol, 33, 865 (1971).
3. Young, D.S., Pestaner, L.C., and Gilberman, U., "Effects of Drugs on clinical Laboratory Tests." Clinical Chemistry, 21, 3660 (1975).
4. Sterling, L., Diagnosis and Treatment of Thyroid Disease, Cleveland CRC press, P.19-51 (1975).

Ed. 01/2015

DCM045-192-4

**DiaMetra S.r.l. Headquarter:** Via Calabria 15

20090 SEGRATE (MI) Italy

Tel. +39-02-2139184

Fax +39-02-2133354

**Manufactory:** Via Pozzuolo 14, 06038 SPELLO (PG) Italy

Tel. +39-0742-24851

Fax +39-0742-316197

E-mail: [info@diametra.com](mailto:info@diametra.com)

DCM045-192-4  
Ed. 01/2015

## T4 ELISA

for routine analysis

Direct immunoenzymatic determination of thyroxine (T4) in human serum or plasma.

IVD

LOT

See external label

2°C 8°C

Σ = 192 test

REF DKO045

### INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of thyroxine (T4) concentration in human serum and plasma.

T4 ELISA kit is intended for laboratory use only.

### 1. CLINICAL SIGNIFICANCE

The thyroid hormone, thyroxine (T4) is produced by the thyroid gland. An important component in the synthesis is iodine. The major form of thyroid hormone in the blood is thyroxine (T4). Thyroxine is converted to the active T<sub>3</sub> (three to four times more potent than T4) within cells by deiodinases (5'-iodinase).

Thyroxine-binding globulin (TGB) is the major carrier protein for circulating thyroid hormone. Only a very small fraction of the circulating hormone is free (unbound) - T4 0.03%.

The thyronines act on the body to increase the basal metabolic rate, affect protein synthesis and increase the body's sensitivity to catecholamines (such as adrenaline) by permissiveness. The thyroid hormones are essential to proper development and differentiation of all cells of the human body. These hormones also regulate protein, fat, and carbohydrate metabolism, affecting how human cells use energetic compounds. Numerous physiological and pathological stimuli influence thyroid hormone synthesis.

Thyroid hormones act through an unknown mechanism to inhibit neuronal activity; one of the effects is the reduction of the body temperature.


Thyrotoxicosis or hyperthyroidism is the clinical syndrome caused by an excess of circulating free thyroxine, free triiodothyronine, or both.

Both T3 and T4 are used to treat thyroid hormone deficiency (hypothyroidism).

### 2. PRINCIPLE

The T4 (antigen) in the sample competes with the antigenic T4 conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti T4 coated on the microplate (solid phase) (the enzyme conjugate should have no measurable binding to serum proteins especially TBG and albumin).

The interaction is illustrated by the following equation:

Ab<sub>cw</sub>: monospecific immobilised antibody (constant quantity)

Ag: native antigen (variable quantity)

EnzAg: antigen conjugated to enzyme HRP (constant quantity)

Ag Ab<sub>cw</sub>: antigen-antibody complex

EnzAg Ab<sub>cw</sub>: antigen-HRP-antibody complex

K<sub>a</sub>: rate constant of association

K<sub>-a</sub>: rate constant of disassociation

K = k<sub>a</sub> / k<sub>-a</sub>: equilibrium constant

After incubation, the bound/free separation is performed by a simple solid-phase washing.

Then, the enzyme HRP in the bound-fraction reacts with the Substrate (H<sub>2</sub>O<sub>2</sub>) and the TMB Substrate and develops a blu color that changes into yellow when the Stop Solution (H<sub>2</sub>SO<sub>4</sub>) is added.

The colour intensity is inversely proportional to the T4 concentration of in the sample. T4 concentration in the sample is calculated through a calibration curve.

### 3. REAGENTS, MATERIALS AND INSTRUMENTATION

#### 3.1. Reagents and materials supplied in the kit

##### 1. Calibrators (6 vials, 2 mL each)

CAL0

REF DCE002/4506-0

CAL1

REF DCE002/4507-0

CAL2

REF DCE002/4508-0

CAL3

REF DCE002/4509-0

CAL4

REF DCE002/4510-0

CAL5

REF DCE002/4511-0

##### 2. Control (1 vial, 2 mL)

Concentration of Control is indicated on Certificate of Analysis

REF DCE045/4503-0

##### 3. Conjugate Buffer (1 vial, 25 mL)

Buffered solution with preservative and binding protein inhibitors

REF DCE002/4501-0

##### 4. Conjugate (1 vial, 2,8 mL)

T4 conjugated with horseradish peroxidase (HRP)

(keep away from sunlight) REF DCE002/4502-0

5. Coated Microplate (2 microplate breakable)  
Antibody anti-T4 adsorbed on microplates  
**REF DCE002/4503-0**
6. TMB Substrate (1 vial, 22 mL)  
H<sub>2</sub>O<sub>2</sub>-TMB 0.26 g/L (avoid any skin contact)  
**REF DCE004-0**
7. Stop Solution (1 vial, 22 mL)  
Sulphuric acid 0.15 mol/L (avoid any skin contact)  
**REF DCE005-0**
8. 50X Conc. Wash Solution (1 vial, 40 mL)  
NaCl 45 g/L; Tween-20 55 g/L **REF DCE006-0**

### 3.2. Reagents necessary not supplied

Distilled water.

### 3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

### Note

Store the reagents at 2-8°C in the dark. Open the bag of reagent 5 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, the microplate is stable until the expiry date of kit.

## 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Proclin 300<sup>R</sup> as preservative. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.
- The concentration of total thyroxine in serum depends on several factors: the function of the thyroid gland and its regulation, thyroxine binding globulin (TBG) and thyroxine binding to TBG. However, the concentration of total thyroxine alone is not sufficient to monitor the clinical status.
- Total serum thyroxine values may increase during pregnancy or administration of oral contraceptives. The table of drugs and interfering conditions in which the values of total thyroxine are affected have been compiled by the Journal of the American Association of Clinical Chemists.
- Not intended for newborn screening

## 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems, it is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

## 6. PROCEDURE

### 6.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>5</sub>)

The serum are ready to use, are calibrated against Human Serum Reference for thyroxine and have approximate concentrations of:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
µg/dL	0	2.0	5.0	10.0	15.0	25.0

**Calibrators concentration is lot-specific; the exact values are stated on labels and Certificate of Analysis for each lot.** Once opened, the Calibrators are stable 6 months at 2-8°C.

### 6.2. Preparation of diluted Conjugate

Dilute the T4 conjugate 1:11 with Total T4 Conjugate Buffer in a suitable container. For example, dilute 160 µL of conjugate in 1.6 mL of conjugate buffer. This reagent should be used within 24 hours for maximum performance of the assay. Stored at 2-8°C.

**IMPORTANT NOTE: prolonged exposure of the Conjugate to sunlight may affect the functional characteristics of the assay; therefore do not expose Conjugate (and diluted Conjugate) to direct sunlight.**

### 6.3. Preparation of Wash Solution

Dilute the contents of each vial of the "50X Conc. Wash Solution" with distilled water to a final volume of 2000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

### 6.4. Preparation of the Sample

The determination of T4 should be performed in human serum or plasma.

Specimens may be refrigerated at 2-8°C (for a maximum period of 48 hours). If the specimens cannot be assayed within 48 hours, the samples may be stored at temperatures of -20°C for up to 30 days. Avoid repetitive freezing and thawing of samples. When assayed in duplicate, 0.050 mL of the specimen is required.

### 6.5. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>), two for each Control, two for each sample, one for Blank.

absorbent paper towel.			
<b>Automatic washer:</b> if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate 15 minutes in the dark at room temperature (22-28°C).			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

## 7. QUALITY CONTROL

Each laboratory should assay controls at levels in the hypothyroid, euthyroid and hyperthyroid range for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

## 8. RESULTS

### 8.1. Mean Absorbance

Calculate the mean of the absorbance (E<sub>m</sub>) for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>) and of each sample.

### 8.2. Calibration curve

Plot the mean value of absorbance (E<sub>m</sub>) of the Calibrators (C<sub>0</sub>-C<sub>5</sub>) against concentration. Draw the best-fit curve through the plotted points. (es: Four Parameter Logistic).

### 8.3. Calculation of Results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in µg/dL.

Reagent	Calibrator	Sample/Control	Blank
Sample/Control		25 µL	
Calibrator C <sub>0</sub> -C <sub>5</sub>	25 µL		
Diluted Conjugate	100 µL	100 µL	
Incubate 1 h at room temperature (22-28°C). Remove the contents from each well, wash the wells three times with 300 µL of diluted wash solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an			

## 9. REFERENCE VALUES

A study of euthyroid adult population was undertaken to determine expected values for the T4 EIA Test System.

	Mean (µg/dL)	SD	Range (µg/dL)
Normal Values	7,6	1,6	4,4 - 10,8

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

## 10. PERFORMANCE AND CHARACTERISTICS

### 10.1. Precision

#### 10.1.1. Intra Assay Variation

Within run variation was determined by replicate measurements (20x) of three different control sera in one assay. The within assay variability is  $\leq 8.16\%$ .

#### 10.1.2. Inter Assay Variation

Between run variations was determined by replicate measurements (10x) of three different control sera in different lots of kit. The between assay variability is  $\leq 8.42\%$ .

### 10.2. Accuracy

The recovery of 1.25 - 2.50 - 5.00 - 10.00 µg/dL of T4 added to sample gave an average value ( $\pm$ SD) of  $97.8\% \pm 5\%$  with reference to the original concentrations.

### 10.3. Sensitivity

The lowest detectable concentration of T4 that can be distinguished from the Calibrator 0 is 0.4 µg/dL at the 95% confidence limit.

### 10.4. Specificity

The cross-reactivity of the triiodothyronine antibody to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated by deriving a ratio between doses of interfering substance to dose of Triiodothyronine needed to displace the same amount of tracer.

Substance	Cross Reactivity	Concentration
I-Thyroxine	1.0000	-
d-Thyroxine	0.9800	10 µg/dL
I-Triiodo-thyronine	0.0300	100 µg/dL
d-Triiodo-thyronine	0.0150	100 µg/dL
Monoiodo-Tyrosine	N/D	100 µg/mL
Diiodo-Tyrosine	N/D	100 µg/mL
Triiodothyroacetic Acid	N/D	100 µg/mL
Tetraiodothyroacetic Acid	0.0001	100 µg/mL

### 10.5. Correlation with RIA

Diametra T4 ELISA was compared to another commercially available T4 assay. 114 serum samples were analysed according in both test systems.

The linear regression curve was calculated:

$$(T4 \text{ Diametra}) = 0.959 \cdot (T4 \text{ RIA}) + 0.339$$

$$r^2 = 0.936$$

## 11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

## BIBLIOGRAPHY

1. Barker, S.B., "Determination of Protein Bound Iodine" *Journal Biological Chemistry*, **173**, 175 (1948).
2. Chopra, I.J., Solomon, D.H., and Ho, R.S., "A Radioimmunoassay of Thyroxine," *J.Clinical Endocrinol.*, **33**, 865 (1971).
3. Young, D.S., Pestaner, L.C., and Gilberman, U., "Effects of Drugs on clinical Laboratory Tests." *Clinical Chemistry*, **21**, 3660 (1975).
4. Sterling, L., *Diagnosis and Treatment of Thyroid Disease*, Cleveland CRC press, P.19-51 (1975).

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**ERROR POSSIBLE CAUSES / SUGGESTIONS****No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

**Too low reaction (too low ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

**Too high reaction (too high ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

**Unexplainable outliers**

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation